




AlModification | No. | Association | No. |





1







HKU LKS Faculty of Medicine
Department of Orthopaer
& Traumatology Current technologies and industries Current technologies 1) Telemedicine > e-medical records, remote care, e-appointment booking, self-symptom checkers, patient outcome reporting, etc. > Teledoc (1st), MeMD, iCliniq, etc. (mostly internal medicine) 2) Wearable devices > smartwatches and on-body sensors for detecting health-related data etc. > AliveCor (mobile ECG), BioTelemetry (health management in clinical setting), Gentag (diagnostic sensors) 3) Augmented reality > technology enhances real-world experiences with computerised sensory information and is used to build smart devices for healthcare professionals > FundamentalVR (surgical training), Karuna Labs (chronic pain management), OxfordVR (Mental Health)

> Assistive technologies, rehabilitation robotics, unobtrusive monitoring sensors, etc.

4) Others





Analysis: Current digital alignment assessment tools Manual alignment analysis product Automated alignment analysis tools Surgimap, X-Align, Integrated Global Alignment, etc. (no product available) Conventional imaging process-based methods and disadvantages · Heterogeneous patterns of deformities X-rays having high variance Limited accuracy<sup>1,2</sup> with original high-resolution images Al integrated methods and disadvantages Semi-automated<sup>3</sup> Automated (CV convenient but not mimic clinicians) · Direct regress CA4 cannot guarantee what has been learnt Indirect regress CA<sup>5</sup> still without intermediate supervision CNN alignment detections<sup>6,7</sup> All used original high-resolution X-rays but often generated low accuracies (>10° errors)

## Our aim and objectives

## Aim


• Provide reliable and easily accessible automated vertebral landmark detection irrespective of image quality.

## Objectives

- Establish a reliable deep learning-based method to accurately detect vertebral landmarks, including endplates and end vertebrae
- Eliminate previous restrictions of automatic coronal alignment on curve patterns or imaging quality by training the model using non-original X-rays of various image quality and different curve patterns.


(G)

10




11 12


\_



Conferent Reference & Try Long to Describe of Odd 10th Description
HEW LASS REQUISITION (Modification Department of Orthopaedics 表 Traumatology 曹港大學城市及顧鄉外科學系


13






15


1



AlignPro™ Automated Alignment: end vertebrae detections Figure 1. Examples of end vertebrae detection Table 1. Evaluation metrics on CA prediction 0.88 0.88 0.88 0.83 0.80 0.82 3.73° 3.64° CA total 0.82 0.83 0.83 4.15° 3.11° · A: CAC and CAT B: CAT and CAL angle C: normal with no curves; . D: CAC and CAT and CAL angle. For each part . (1) the ground truth heatmap · (2) predicted heatmap • (3) original image merged with predicted heatmap An interesting false negative in the cervicothoracic region was shown in D. 1


17





19

\_



Take home messages Clinical relevance key points 1) Al powered Digital Health will be the future of healthcare 2) AlignPro™ using our SpineHRNet enables fully automated alignment analysis with no limitations of the curve patterns using smartphones, providing consistent and fast alignment analysis results for large clinical trials and may facilitate out of hospital consultation. Contact
Nami: Eng Grace Zhang; Jason Cheung | Email: tgzhang@hku.hk; cheungjp@hku.hk |
Phone: 391-76989 | Fax: 281-85210 | Location: level 5 Professorial Building, Queen Mary Hospital

22